


# RP271

## X-Ray Protection for Paediatrics



Size Adjustment From Behind

Article Code: **RP271** Mammary Gland Shield for Children

CE

Made in Germany


### RP271 – X-Ray Protection for Paediatrics

The primary function of this paediatric X-ray protective cover is to protect the mammary gland and surrounding tissue during X-ray images in medical diagnostics.

When correctly placed on the patient, the RP271 can reduce the radiation exposure up to 80%.

The size of this reusable product can be easily adjusted with the comfortable, soft neoprene stoppers in the front and a convenient cord stopper in the back.

- **Protective Material** NovaLite radiation protective material
- **Lead Equivalent** Pb 0.50 mm
- **Outer Cover/Colour** ComforTex® HPMF available in Teddy-Design
- **Size** Universal, easily adjustable



Width approx.  
10 cm

#### MAVIG GmbH

X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone: +49 (0) 89 / 420 96-0  
Fax: +49 (0) 89 / 420 96-200  
e-Mail: info@mavig.com

[www.mavig.com](http://www.mavig.com)