

RP643

MAVIG



Article Code: **RP643** Gonadal Apron



AW417 - Wall Mount

CE

Made in Germany

Gonadal Apron RP643

The number of patients passing through today's X-ray facilities is very high making it important for the X-ray shields to be practical, and easy to put on and take off quickly. The flexible tension strap of our series RP643 offers the easiest solution.

Our gonadal aprons are made of high quality, flexible lead rubber multi-ply material with outstanding weight distribution in this material category.

We recommend our specially designed wall rack model AW417 for proper storage of your gonadal aprons.

- **Lead Equivalent** Pb 0.50 mm - multilayer lead rubber material
- **Sizes**
 - S - for small children, covers 20 x 25 cm (HxW)
 - M - for children, covers 30 x 30 cm (HxW)
 - L - for adults, covers 37 x 40 cm (HxW)
 - XL - for adults, covers 45 x 50 cm (HxW)
- **Outer Cover/Colour** Comfortex® HPMF, available in the colours:
Curacao, Limette, Regatta, Indian Summer, Ocean, Lagoon, Orchidee, and Teddy design

MAVIG GmbH

X-Ray Protection
& System Solutions

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone: +49 (0) 89 / 420 96-0
Fax: +49 (0) 89 / 420 96-200
e-Mail: info@mavig.com

www.mavig.com