

WD255




MAVIG

Mobile X-Ray Protection WD255

Solid X-Ray protection can also be economical.

That is why the WD255 was created as a price sensitive solution ideal for workplace shielding. The slim and light-weight design protects the user from radiation while allowing the shield to take up a minimal amount of space and making it easy to move when not in use.

MAVIG offers the mobile shield WD255 in two different lead equivalent values depending on the client's needs. Customized solutions can also be provided.



Standard Shield Variations, Window (30 x 40 cm)

Article No.	Size (W x H)	Lead Value
WD2551B	100 x 200 cm, with window	Pb 1.00 mm
WD2552B	100 x 200 cm, with window	Pb 2.00 mm

