


EN

# BR118

**Designs** Taupe  
Blue  
Black

MAVIG



Taupe



Blue



Black

**Model:** BR118

Glasses with front X-ray protection of Pb 0.75 mm, universal size, weight approx. 60 g

Made in Germany



## BR118 – X-Ray Protective Glasses

The X-ray protective glasses BR118 have curved glasses and a front protection of Pb 0.75 mm. The frame is made of an extremely durable plastic that ensures easy cleaning and good resistance to chemicals (disinfectant agents). The BR118 comes with a storage case and a soft google strap.

MAVIG uses only protective lenses of the highest quality. This guaranties the best light transmission values and a maximum lens break resistance. An ideal composition of the protective material results in a high level X-ray absorption, and optimized production processes achieve an absolute homogeneous allocation of the glass additives.

The quality features of our glasses make them perfect for the use in the medical environment:

- Great frame stability through high-grade synthetic materials
- Optimized balance of sturdy, light-weight frame and lenses
- Minimized risk of injury due to shape and used components
- Low allergy risk due to carefully chosen materials
- Easy cleaning and superior resistance against chemical influences (disinfection)
- For this model prescription lenses are available: single and bifocal

**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail [info@mavig.com](mailto:info@mavig.com)