


APRON FOR FRONT PROTECTION



COMFORT RA660

This model complies with IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016.

Naturally, MAVI G X-ray protective clothing also underwent the required prototype testing and certification in compliance with PPE 89/686/EEC (CE0123).

Apron for Front Protection

COMFORT RA660

Comfort RA660

Up-to-date radiation protective clothing which accompanies your active working day and provides outstanding qualifications.

Our customers declare that this model is their favourite front apron. It stands out because it can be put on and taken off easily, fits the body well, stays in place and provides outstanding shoulder relief.

The exceptional comfort of the Comfort RA660 series is enhanced by the wide stretch inserts on the back wings.

MAVIG X-ray protective materials achieve an outstanding rating in terms of weight and X-ray protection. Produced with the first-class X-ray protective material NovaLite, the apron ensures fulfilment of fundamental occupational health requirements for long wearing periods in radiological workplaces.

The RA660 model is now equipped with shoulder pads as standard for optimised wearing comfort.



Special Features

The wide stretch inserts in the back wings ensures individual adjustment, a good fit and shoulder relief, even when the X-ray protective clothing is worn for a long time.

To reduce the weight of the apron, "without protection" can be optionally chosen in the back area.

Closure System

MAVIG uses only top quality hook & burr fasteners.



If wear and tear should occur on the fleece of our hook & burr closure, it can easily be exchanged.

Protective Material

- NovaLite or
- Standard Lead

Lead Equivalent Values

Front side, optionally: Pb 0.25 mm,
Pb 0.35 mm, or
Pb 0.50 mm
Back, optionally: Pb 0.25 mm or
without protection

Outer Material / Colours

Optionally ComforTex® HPMF Hybrid or ComforTex® HPMF, available in the colours:

- Curacao
- Orchidee
- Regatta
- Indian Summer
- Limette
- Ocean
- Lagoon

Sizes*, recommended for:

- XS** Chest girth ≤ 84 cm
Waist/hip girth ≤ 91 cm
- S** Chest girth approx. 85 - 95 cm
Waist/hip girth approx. 92 - 102 cm
- M** Chest girth approx. 96 - 106 cm
Waist/hip girth approx. 103 - 112 cm
- L** Chest girth approx. 107 - 117 cm
Waist/hip girth approx. 113 - 123 cm
- XL** Chest girth approx. 118 - 128 cm
Waist/hip girth approx. 124 - 134 cm
- XXL** Chest girth approx. 129 - 139 cm
Waist/hip girth approx. 135 - 145 cm
- 3XL** Chest girth approx. 140 - 150 cm
Waist/hip girth approx. 146 - 156 cm

* The size depends on the girth of the chest, waist, or hip, depending on which of the three measurements is the largest.

Lengths, recommended height:

- 90 cm for approx. 150 - 160 cm
100 cm for approx. 160 - 170 cm
110 cm for approx. 170 - 180 cm
120 cm for approx. 175 - 190 cm
130 cm for approx. 185 - 195 cm
140 cm for approx. 195 - 205 cm

Optional / On Request

Sleeves

Breast Pocket

Embroidery:

- on the breast pocket
- on permanently sewn-on border
- on detachable border

Thyroid Protection

Made-to-Measure Protection
(possible on an individual basis)

Replacement Order: Elastic Bands with a Hook & Burr Closure (fleece)
(adapted via zipper)