



# EN RD614

MAVIG



**Model Variants:** RD614E  
RD614K  
RD614EAK  
RD614KAK

Patient Thyroid and Sternum Protection with hook & burr closure for adults  
Patient Thyroid and Sternum Protection with hook & burr closure for children  
Similar to RD614E with an exchangeable hook & burr closure for adults  
Similar to RD614K with an exchangeable hook & burr closure for children

Made in Germany



## RD614 – Patient Thyroid and Sternum Protection

The thyroid is a radiation-sensitive organ of the human body which requires particular protection.

Our thyroid and sternum protector is characterised by its first-class cut and the best fit, taking into account all aspects of radiation protection.

A simple, optionally replaceable hook & burr fastener ensures sufficient room for adjustment.

- **Protective Material** Standard lead
- **Lead Equivalent** Pb 0.50 mm
- **Outer Material / Colours** ComforTex® HPMF, available in the colours:  
Curacao, Limette, Orchidee, Ocean, Regatta, Lagoon, Indian Summer, Teddy design
- **Neck Sizes** Standard size fits most neck sizes,  
RD614E: 32 - 45 cm / RD614EAK: 37 - 50 cm / RD614K & RD614KAK: 28 - 36 cm
- **Hygienic Cover** Washable hygienic cover made of the special ComforTex® MF material in the colour "Ocean"  
(article codes: RD614EF - adults,  
RD614KF - children)



**MAVIG GmbH**  
X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone +49 (0) 89 / 420 96 - 0  
Fax +49 (0) 89 / 420 96 - 200  
e-Mail info@mavig.com