

# RD635



## Dental Apron for Images when Standing

- **Lead Equivalent:**  
Pb 0.50 mm
- **Inner Material:**  
Multi-Layer Lead Rubber
- **Outer Cover:**  
ComforTex® HPMF
- **Size:**  
Children (K), Adults (E)

**RD635ECC\*** X-Ray Protective Dental Apron for Adults  
**RD635KCC\*** X-Ray Protective Dental Apron for Children

Dimensions: Length 90.0 cm; Width 60.0 cm  
Dimensions: Length 70.0 cm; Width 45.0 cm

MAVIG



CE

Made in Germany

## MAVIG X-Ray Protective Dental Apron

The patient protective apron RD635 offers ideal radiation protection during panoramic X-ray imaging, cephalometry, or cone beam computed tomography (CBCT) of a standing patient.

The apron can be put on comfortably, and with only one movement, it can be individually adjusted to the specific measurements of the respective patient.

Another advantage is the hook & burr closing system in the front of the apron, which allows the back part of the apron to lay flat against the patient.

### Inner Material: Multi-Layer Lead Rubber with Pb 0.50 mm

Our natural lead rubber, which does not contain any artificial plasticisers, has kept us in the lead as a manufacturer of X-ray protective clothing for many years, with excellent marks for weight distribution and radiation protection.

### Outer Cover: ComforTex® HPMF - Single and Dual Colour Designs

Our ComforTex® HPMF, a high performance medical fabric according to DIN EN 13795, provides a far greater protection of the radiation protective inlay than ever before. Its improved barrier function complies with the very demanding hygienic standards in today's work environment.

\*CC - ColourCode: 1 = Curacao, 2 = Orchidee, 3 = Regatta, 4 = Indian Summer, 5 = Limette, 6 = Teddy Design, 7 = Ocean, 8 = Lagoon

### MAVIG GmbH

X-Ray Protection and  
Medical Suspension Systems

Stahlgruberring 5  
81829 Munich  
Germany

PO Box 82 03 62  
81803 Munich  
Germany

Phone: +49 (0) 89 420 96 - 0  
Fax: +49 (0) 89 420 96 - 200  
e-Mail: info@mavig.com

[www.mavig.com](http://www.mavig.com)