

RD642

MAVIG

Dental Apron

for Images
when Sitting
or Lying

- **Lead Equivalent:**
Pb 0.50 mm
- **Inner Material:**
Multi-Layer Lead Rubber
- **Outer Cover:**
ComforTex® HPMF
- **Size:**
Children (K), Adults (E)



RD642ECC* X-Ray Protective Dental Apron for Adults
RD642KCC* X-Ray Protective Dental Apron for Children

Dimensions: Length 90.0 cm; Width 60.0 cm
Dimensions: Length 70.0 cm; Width 45.0 cm

CE

Made in Germany

MAVIG X-Ray Protective Dental Apron

The MAVIG patient protective apron RD642, which is used for dental images in both sitting and lying positions, distinguishes itself with its flexibility and wearing comfort.

High quality hook and burr closures enable precise fitting, as well as safe and fast putting on of the apron.

The apron extensively protects the front of the sitting or lying patient against radiation.

Inner Material: Multi-Layer Lead Rubber with Pb 0.50 mm

Our natural lead rubber, which does not contain any artificial plasticisers, has kept us in the lead as a manufacturer of X-ray protective clothing for many years, with excellent marks for weight distribution and radiation protection.

Outer Cover: ComforTex® HPMF - Single and Dual Colour Designs

Our ComforTex® HPMF, a high performance medical fabric according to DIN EN 13795, provides a far greater protection of the radiation protective inlay than ever before. Its improved barrier function complies with the very demanding hygienic standards in today's work environment.

MAVIG GmbH

X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone: +49 (0) 89 420 96 - 0
Fax: +49 (0) 89 420 96 - 200
e-Mail: info@mavig.com

*CC - ColourCode: 1 = Curacao, 2 = Orchidee, 3 = Regatta, 4 = Indian Summer, 5 = Limette, 6 = Teddy Design, 7 = Ocean, 8 = Lagoon

www.mavig.com