

RP270 Breast Protection



Article Code: **RP270** Breast Protection for Women


CE
Made in Germany

RP270 Breast Protection for Women

MAVIG's breast protection was specially designed for the protection of the mammary gland area of the patients during mammographic radiation procedures. As all of our protection the RP270 offers solid protection and optimized functionality.

The protector consists of two parts, each 18 x 40 cm (W x H), joined by an adaptable hook & burr closure for a comfortable fit and easy size adjustment.

- **Lead Equivalent** Pb 0.50 mm - NovaLite X-ray protective material
- **Size** Universal
- **Outer Cover/Colour** ComforTex® HPMF, available in the colours: *Curacao, Limette, Regatta, Indian Summer, Ocean, Lagoon, Orchidee and Teddy Design*
- **Edging** All edges are covered with a seem band.



MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone: +49 (0) 89 / 420 96-0
Fax: +49 (0) 89 / 420 96-200
e-Mail: info@mavig.com

www.mavig.com