



RP278



MAVIG

Placement: Spina iliaca ventralis



Article Numbers:	RP2781	Ovarian Protection, Large
	RP2782	Ovarian Protection, Medium
	RP2783	Ovarian Protection, Small
	RP278SET	Ovarian Protection Set, all three sizes

CE

Made in Germany

RP278 - Ovarian X-Ray Protection for Female Patients

The development of the MAVIG ovarian protection is unrivalled and meets all requirements of a radiology practice. The radiation exposure to the gonads of female lying patients, for example during pelvic exams, is significantly lowered without hindering the diagnostic technique.

The two radiolucent control marks ensure correct and uncomplicated placement. In addition to the three available sizes, each shield can be positioned by means of two swiveling parts to allow for fine adjustment.

Of course, the RP278 is easy to clean and disinfect.

- **Outer Material** Soft PVC Cover (light gray colour)
- **Protective Material** Multi-Layer Lead Rubber
- **Lead Equivalency** Pb 1.00 mm
- **Sizes**
 - Large – Height 12 cm
 - Medium – Height 9 cm
 - Small – Height 7 cm

* Available individually or in a set of all three sizes.

MAVIG GmbH
X-Ray Protection
& System Solutions

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone: +49 (0) 89 / 420 96-0
Fax: +49 (0) 89 / 420 96-200
e-Mail: info@mavig.com

www.mavig.com