

EN

RP679

MAVIG



Female Patient Ovaries Protection with Belt

Sets:	RP679SET-E	Ovarian protection with stretch belt RP679 for adults (set with 3 sizes)
	RP679SET-K	Ovarian protection with stretch belt RP679 for children (set with 3 sizes)

Made in Germany



RP679 – Ovarian Protection with Stretch Belt

The patient protection RP679 protects the ovaries of your female standing patients during X-rays in medical diagnostics (set RP679SET-E). The elastic stretch belt with snap buckle guarantees easy putting on and taking off and a comfortable size adjustment.

With the set RP679SET-K the ovarian protection RP679 is also available as patient protection for paediatrics.

- **X-ray protective material Standard-Lead with a lead equivalent of Pb 1.00 mm**
- **Comfortable size adjustment by means of the stretch belt with hygienic cover and snap buckle**
- **Quick and easy to put on**
- **Set consisting of the stretch belt* and three protectors with the following dimensions:**
for adults: Small (Ocean, 16 x 11 cm), Medium (Indian Summer, 18 x 13 cm), Large (Regatta, 20 x 19 cm)
for children: Small (Orchidee, 7 x 6 cm), Medium (Limette, 9 x 8 cm), Large (Curacao, 10 x 10 cm)
- **Protectors can be attached to the belt with the hook & burr fastening**

* A custom size belt is available for patients with a girth greater than 128 cm.

MAVIG GmbH
X-Ray Protection and
Medical Suspension Systems

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone +49 (0) 89 / 420 96 - 0
Fax +49 (0) 89 / 420 96 - 200
e-Mail info@mavig.com