



PORTEGRA²

Wall Mounting

TS1515 Wall mount for Portegra2
with a max. load capacity of 13.0 kg (28.6 lbs)


CE, IEC 60601-1
Made in Germany

MAVIG Wall Mounting for Portegra2

MAVIG has set a standard for practices and hospitals in ceiling support and radiation protective systems.

In situations where installing a ceiling column is not an option, MAVIG offers a second possibility with the wall mount for Portegra2. Just like the ceiling suspension systems, the wall mount allows the user to conveniently free up floor space for other medical devices without large investment.

- **One mounting pin**
- **Easy installation**
- **Cable cut-out for external cabling of electrical equipment**
(pin is not electrified)
- **Max. load capacity at spring arm adapter:**
13.0 kg (28.6 lbs) for 750 mm extension arm
10.0 kg (22.0 lbs) for 950 mm extension arm



MAVIG GmbH
X-Ray Protection
& System Solutions

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone: +49 (0) 89 420 96 - 0
Fax: +49 (0) 89 420 96 - 200
e-Mail: info@mavig.com

www.mavig.com