

Twin

360°

PORTEGRA²

Ceiling Columns for Track Mounting

TS2031 - 32 360° column with carriage (please refer to table below for different lengths)
TS1501 - 03 Twin column with carriage (please refer to table below for different lengths)



CE, IEC 60601-1
 Made in Germany

MAVIG Columns with Carriage – 360° and Twin

Solid construction for maximum safety. Many years of development has led to MAVIG's unique steel columns. The proven construction, sturdy design, and numerous tests ensure the highest level of safety available. The patented brake mechanism is a standard for all support systems with a carriage. The brake holds the carriage at a fixed position inside the ceiling track, which provides additional safety during its use.

MAVIG columns are equipped standard with one electrical pin for a direct connection to a lamp.

- Each pin offers a 240° rotation (360° for the lower pin of the 360° column)
- One electrical connection: One 3-pole sliding contact
- Maximum net load at spring arm adapter: 18.0 kg (39.7 lbs) per pin
- Different lengths to match room height

MAVIG GmbH
 X-Ray Protection
 & System Solutions

Stahlgruberring 5
 81829 Munich
 Germany

PO Box 82 03 62
 81803 Munich
 Germany

Phone: +49 (0) 89 420 96 - 0
 Fax: +49 (0) 89 420 96 - 200
 e-Mail: info@mavig.com

Article-No. 360°	Length*	Article-No. Twin	Length*
TS2031	800 mm	TS1501	800 mm
TS2032	570 mm	TS1502	580 mm
		TS1503	460 mm

* Ceiling column lengths are measured from the top of the mounting flange to the bottom of the column structure. For a total length when installed, 75 mm needs to be added for the track and 60/80 mm for the brake handle.

