


PORTEGRA²

TriPin Ceiling Column

*Stationary
Mounting*



TS2320 - 22 TriPin column stationary mounting with safety steel rope (see table below for lengths)



Made in Germany

MAVIG Portegra2 - TriPin Column

The innovative TriPin ceiling column series expands on Portegra2's successful slim design and flexibility to provide a cost-effective way to manage three devices within a compact space. Portegra2's signature 360° mounting pin still guarantees the highest degree of moveability and positioning.

MAVIG TriPin columns are equipped standard with two electrical connections for fixed mounted systems.

- **Three mounting posts**
- **Both upper pins offer a 240° rotation; Lower pin 360°**
- **Electrical connection: 3-pole sliding contacts**
- **Maximum 18.0 kg (39.7 lbs) net load for each mounting post**
- **Different lengths to match room height**
- **Compatible with existing Portegra2 mounting plates and suspension arms**

Article-No. (stationary)	Length*
TS2320	1000 mm
TS2321	850 mm
TS2322	580 mm

* Ceiling column lengths are measured from the top of the mounting flange to the bottom of the column structure.

MAVIG GmbH
X-Ray Protection
& System Solutions

Stahlgruberring 5
81829 Munich
Germany

PO Box 82 03 62
81803 Munich
Germany

Phone: +49 (0) 89 420 96 - 0
Fax: +49 (0) 89 420 96 - 200
e-Mail: info@mavig.com

www.mavig.com