



EN

WD261

MAVIG



Trouble-Free Movement

Article Number: **WD26105**
WD26105-S
 STEA-WD261

Up close, height adjustable mobile shield, incl. 1 box of sterile covers
 Short version: Up close, height adjustable mobile shield, incl. 1 box of sterile covers
 Box of 20 sterile covers for WD26105 or WD26105-S

Made in Germany



WD261 – Radiation Protective Mobile Shield

Radiation protection that does not affect your safety and work routine

The compact, anatomically formed shield of the model WD261 allows for a close body fit so that the doctor can move freely throughout the procedure. Due to springs the user can adjust the height of the shield to his own height or working position requested, at any time i.e. also during the procedure.

The shield was designed to offer sufficient legroom in order to eliminate any risk of injury. The unimpeded use of the control pedal of the x-ray unit accounts for another advantage.

- **Height adjustable for individual needs from 1060 to 1360 mm** (short version: 910 to 1210 mm)
- **Lead equivalent Pb 0.50 mm**
- **Easy to clean outer material: MAVIG ComforTex® HMPF**
- **Designed for easy and comfortable control pedal use**
- **High quality, electrically conductive wheels with brakes**
- **Specially designed sterile covers STEA-WD261 to meet strict hospital hygienic regulations**

MAVIG GmbH

X-Ray Protection and
 Medical Suspension Systems

Stahlgruberring 5
 81829 Munich
 Germany

PO Box 82 03 62
 81803 Munich
 Germany

Phone +49 (0) 89 / 420 96 - 0
 Fax +49 (0) 89 / 420 96 - 200
 e-Mail info@mavig.com