

STEA-RZAM – Cover for Radiation Protective Drape Radial

Sterile Covers suitable for the radiation protective drape for radial access.

Characteristics

The sterile covers allow for a positioning of the drapes above the sterile patient cover and thus the easy removal in case of an emergency. The special material developed for these covers assures that the drape will not slide during the intervention performed.


Technical Data

Scope of Delivery: Box with 50 Covers

Dimensions 41 x 86 cm (W x H)

Compatible with:

- Radiation Protective Drape ST-RZ5AMM / ST-RZ5AMT



Are you interested in our steril cover STEA-RZAM? Please send us a nonbinding request for a quote.

ADD TO CART