

UT6120 – Lower Body Protection Right of Adapter

The UT61 is MAVIG's answer to table mounted protection for tiltable tables. The lower sway panels self adjust to a table inclination of 15°. The adapter is on the left side of the lamella.

Characteristics

The lower part consists of three overlapping sway panels, the upper shield of one panel reinforced by a 3 mm PVC board. All panels are made of flexible multilayer lead rubber with a light grey and hygienic PVC covering.

All metal parts (adapter etc.) are made of high quality stainless steel or anodized Aluminium.

The scope of delivery includes a wall mount for practical and space saving storage if the shield is not used.

Technical Data

Scope of Delivery: Lower part with adapter, upper shield, wall mount

Dimensions	Lower Part	50 x 77 cm (W x H)
	Upper Shield	50 x 30 cm (W x H)
	Adapter	suitable for flat profile rails 10 x 25 mm
Weight	Lower Part	6.00 kg
	Upper Shield	2.50 kg
	Wall mount	2.90 kg
Lead Equivalent	Pb 0.50 mm	

Compatible with:

- Artis zee MP



Are you interested in our UT6120? Please send us a nonbinding request for a quote.

ADD TO CART