

MD-69 – Mobile Accessories Cart

The mobile accessories cart MD-69 from MAVIG can easily hold any device which can be attached to a rail; such as the control panels of a C-arm. The cart can be moved to different locations in the room, for example to maximize the distance of the medical staff from the source of x-ray radiation.

Characteristics

The MD-69 has been especially designed to require minimal space while still providing maximum tip resistance. The mobile accessories cart has a 72 cm long 10 x 25 mm flat profile DIN rail with a load capacity of up to 8.5 kg which allows you to easily mount various devices.

For effortless maneuverability, the cart has been equipped with smooth rolling, electroconductive castors. To secure the cart in place, the front two castors are fitted with brakes.

Technical Data

Scope of Delivery: Mobile Accessories Cart

Dimensions (H x B x T)	DIN-Rail (H x B x T)	Max. Load
94. x 96.0 x 64.4 cm	25 x 720 x 10 mm	8.5 kg

Compatible with:

- Ceiling Suspended Angiography System: Artis zee, Artis Q, Artis Q.zen
- Floor Mounted Angiography System: Artis zee, Artis Q, Artis Q.zen
- Biplane Angiography System: Artis zee, Artis Q, Artis Q.zen
- Artis zeego (Ceiling Height > 2.90 m)
- Artis one



Are you interested in our Mobile Accessories Cart MD-69? Please send us a nonbinding request for a quote.

ADD TO CART