

## UT3075SET – Radiation Protective Panels

A solid radiation protection concept should be an integral part of any radiological procedure. The UT5001-SI offers perfect protection especially in the gastrointestinal work environment. In order to extend the lateral protection on the table serves the set UT3075, consisting of three individual panels.

Table mounted shields substantially reduce scattered radiation to the medical personnel, especially to the lower extremities. For increased safety, MAVIG's innovative table mounted, protective shields are solely produced from materials with the highest quality.

### Characteristics

MAVIG's system UT30 makes it possible to attach panels to the table individually. In a single step they are fixed securely and overlapped for safety.

The radiation protection panels synchronize to all movements of the table, effectively protecting the medical staff.

### Technical Data

**Scope of Delivery:** 3x Radiation Protective Panel

Component	Dimension (LxW)	Lead equivalent	Weight
Panel	je 75.0 x 29.0 cm	Pb 0.50 mm	je 1.75 kg / Set

**Compatible with:** Artis zee MP



**Are you interested in our Radiation Protective Panels UT-3075SET? Please send us a nonbinding request for a quote.**

**ADD TO CART**