

UT5001-SI – Extended Lower Body Protection System

A solid radiation protection concept should be an integral part of any radiological procedure. The UT5001-SI offers perfect protection especially in the gastrointestinal work environment. Physicians and medical staff will find a solution optimized for their application which does not interfere with their workflow.

Table mounted shields substantially reduce scattered radiation to the medical personnel, especially to the lower extremities. For increased safety, MAVIG's innovative table mounted, protective shields are solely produced from materials with the highest quality.

Characteristics

The lower shield sections are comprised of overlapping protective sway panels, which self-adjust to different table inclinations of up to 15°. The new multi-functional side handles allow the user to easily carry, position, and firmly lock the lower shields in place.

The various top shields of the UT5001-SI can be combined and strategically placed to allow for multiple interventional procedures. For comfortable patient positioning or in case of emergency, the shields can be removed quickly.

Technical Data

Scope of Delivery: Frame, Lower Front Shield, Lower Side Shield, Top Front Shield, 2x Top Side Shield/Long, 2x Top Side Shield/Short, 2x Wall Mount Lower Shields, 2x Wall Mount Top Shields

Parts	Dimensions (WxH)	Lead Equivalent	Weight
Lower Protection Front	725 x 684 mm	Pb 0.50 mm	6.3 kg
Lower Protection Side	1130 x 684 mm	Pb 0.50 mm	9.8 kg
Top Shield Front	703 x 301 mm	Pb 0.50 mm	3.2 kg
Top Shield Side Long	186 x 351 mm	Pb 0.50 mm	0.8 kg
Top Shield Side Short	186 x 301 mm	Pb 0.50 mm	0.6 kg
Frame	612 x 746 mm	-	7.6 kg

Compatible with: Artis zee MP



Are you interested in our shield OT94001? Please send us a non-binding request for a quote.

[ADD TO CART](#)